
**Final Report**

*Title:* Evaluation of Toxicity Responses of EpiAirway™ to Nutriair®  
Product Exposures

*MatTek Study Number:* 012-18

*Sponsor:* Joshua Matzkin  
NV Nutrition LLC  
3226 Bennett Street North  
St. Petersburg, FL 33713

*Study Director:*

  
Anna Maione, Ph.D.

*Date:*

April 4, 2018

**TABLE OF CONTENTS**

|                                |    |
|--------------------------------|----|
| 1. Quality Assurance Statement | 3  |
| 2. Executive Summary           | 4  |
| 3. Experimental Methods        | 4  |
| 4. Test Articles               | 5  |
| 5. QC Results                  | 6  |
| 6. Experimental Results        | 6  |
| 7. Conclusion                  | 12 |

**Quality Assurance Statement**

I have reviewed the procedures and data presented in this report.

The AIR-100-PC12 tissues (Lot # 28024) and reagents used within this study were produced or obtained in accordance with MatTek's SOPs. The tissue production and test protocols were performed by qualified laboratory personnel.

Based on the results of standardized Quality Control testing, all tissues used in this study fulfill the acceptance criteria for the tissue models.

I hereby certify that the data contained herein are in agreement with the raw data collected during the performance of the study.

*Paul Kearney*

Paul Kearney  
Quality Assurance Director

*April 4, 2018*


Date

## 2. Executive Summary


The current study was conducted to compare the toxicity response of EpiAirway™ tissues to aerosol from Nutriair Energy (20 puffs) and Nutriair Sleep (5 puffs), along with vehicle controls (pure 50:50 PG/VG, Energy vehicle + flavor and Sleep vehicle + flavor). Tissues were exposed in triplicate to aerosol and control air using a Vitrocell VC1 smoking machine. Untreated, incubator control tissues and Triton X-100 treatment control tissues were also included. Aerosol was generated using the CRM81 regime under a 0.5 L/min dilution rate and 20 mL/min vacuum rate. Twenty four hours following exposure, tissue viability was assessed by the LDH assay, barrier function was assessed by TEER and tissue morphology was assessed by H&E staining of tissue cross-sections. All tissues remained viable following exposure to either the full formulation Energy or Sleep aerosol. Likewise, tissue viability was unaffected by exposure to any of the vehicle controls. Similarly, tissue barrier function remained intact after all test article exposures. Histological assessment of tissue sections demonstrated that tissue morphology was unaffected by test article exposure, with all tissues exhibiting a pseudostratified and well-differentiated epithelial layer, stereotypical of the human respiratory epithelium. Taken together, these data indicate that exposure to Nutriair Energy Inhaler or Nutriair Sleep Inhaler did not negatively impact the viability, barrier function or morphology of EpiAirway tissues under the conditions tested. Follow-up experiments could be conducted to confirm and further expand on these results.

## 3. Experimental Methods


The EpiAirway™ model (AIR-100-PC12) is a highly differentiated *in vitro* human airway culture derived from primary human tracheal/bronchial epithelial cells. The human airway epithelial cells have been cultured to form a 3-dimensional model of human airway epithelial tissue. Morphologically, EpiAirway™ is of uniform thickness and is very similar to native nasal and tracheal epithelial tissue in that it exhibits a pseudostratified morphology and contains both ciliated and mucin producing cells. EpiAirway™ tissues possess *in vivo*-like barrier properties due to the formation of functional tight junctions between adjacent epithelial cells. EpiAirway™ cultures also express numerous *in vivo*-like drug metabolizing capabilities including P450, UDP-glucuronosyltransferase, glutathione S-transferase and alkaline phosphatase activities (See MatTek Technical References #456 and 566 and internal data). The EpiAirway™ PC12 tissues (surface 1.12 cm<sup>2</sup>) are cultured on collagen coated microporous membrane cell culture inserts and are grown at the air liquid interface, which allows direct apical application of test materials (similar to *in vivo* exposure). These inserts are the correct size for aerosol exposure using the Vitrocell VC1 smoking machine and 6 well exposure chamber. The model system is shown schematically in **Figure 1**. H&E stained paraffin sections of the EpiAirway™ cultures compared to excised human bronchial epithelium is shown in **Figure 2**. A transmission electron micrograph showing cilia and tight junctions at the apical surface the EpiAirway™ model is shown in **Figure 3**.

**Figure 1.** The EpiAirway Model System showing apical compartment, epithelial tissue, microporous membrane and basal (receiver) compartment.

**Figure 2.** Pseudostratified mucociliary morphology of the EpiAirway tissue compared to excised human bronchial epithelium.

**Figure 3.** Transmission electron micrograph of EpiAirway showing cilia (A) and tight junctions (B).

A detailed protocol is found in **Appendix B**. Briefly, the apical surfaces of EpiAirway™ tissues were rinsed with PBS prior to exposure. Triplicate tissues were exposed to 20 puffs of Energy Inhaler, 20 puffs of Energy Inhaler flavor vehicle control, 5 puffs of Sleep Inhaler, 5 puffs of Sleep Inhaler flavor vehicle control or 20 puffs of PG/VG vehicle control in parallel with exposure to clean control air. All aerosol was generated by a Vitrocell VC1 smoking machine using the CORESTA Recommended Method Number 81(CRM81) regime (square-wave profile, 55 mL puff volume, 3s duration, 8s exhaust and 30s interval) under a 0.5 L/min dilution rate and 20 mL/min vacuum rate. Additional controls, performed in triplicate on each day of Inhaler exposure, included untreated incubator control tissues (negative control) and tissues treated apically for three hours with 0.5% Triton X-100 (positive control for tissue death).

EpiAirway™ tissues were analyzed for barrier function using trans-epithelial electrical resistance (TEER) before and after exposure. Tissue viability was assessed after exposure using the LDH assay and tissue morphology was assessed by Hemotoxylin and Eosin (H&E) staining of tissue sections. In addition, conditioned media were collected from all tissue cultures for possible future cytokine analysis. The mean  $\pm$  standard deviation of TEER values (expressed as percent of pre-exposure value) and percent viability were calculated relative to matched controls. Additional analyses are available upon request.

#### 4. Test Articles

Test articles were supplied by the Sponsor. Triton X-100 was supplied by MatTek. Air used for clean air controls was compressed breathing quality air composed of 76.5-80.5% nitrogen and 19.5-23.5% oxygen.

Additional components are outlined below:

| Article              | CAS #     | Manufacturer/ Catalog # | Lot#                |
|----------------------|-----------|-------------------------|---------------------|
| TritonX-100          | 9002-93-1 | Sigma / X100            | 112217OOB           |
| LDH cytotoxicity kit | N/A       | Takara / MK401          | AH4P029,<br>AH4P028 |

**Table 1. Test Articles.**

Vehicle controls: 50:50 PG/VG, Energy Inhaler vehicle + flavor (Vfrost 50:50 PG/VG), Sleep Inhaler vehicle + flavor (Earl Grey 50:50 PG/VG)

Negative controls: Clean air, Incubator control

Positive controls: 0.5% Triton X-100 (positive control for cytotoxicity)

## 5. Quality Control

Baseline tissue barrier function by TEER:

Study QC acceptance criteria:  $>300 \Omega \cdot \text{cm}^2$

| Lot # | TEER $\pm$ s.d.   |
|-------|-------------------|
| 28024 | 755.9 $\pm$ 153.4 |

**Table 2. QC TEER Data.**

**Conclusion:** Based on the TEER results obtained from the standardized quality control tests, the EpiAirway™ tissues used in this study meet the QC acceptance criteria.


## 6. Results

### 6a. Macroscopic observations of EpiAirway™ tissues following exposure to test articles.

The apical surfaces of tissues exposed to test articles appeared similar to that of the matched control air-exposed tissues and to the untreated, incubator control tissues. Triton-treated tissues were leaky and media had seeped into the apical side of the insert.

### 6b. EpiAirway™ tissue viability determined by LDH following exposure to test articles.

Tissue viability following control treatments and aerosol/air exposure was measured using an LDH assay. Viability was calculated relative to the untreated control (100% viable) and Triton X-100 treatment (kill control). All tissues were  $> 90\%$  viable following treatment with either test articles or control air (**Figure 4a**). There was no difference in viability between the test article-exposed tissues and the matched control air-exposed tissues, indicating that the test articles did not induce any cytotoxicity (**Figure 4b**).


**Figure 4. Tissue viability determined by LDH.** Tissue viability following exposure to test article aerosol or matched control air was assessed using LDH. (a.) Tissue viability calculated relative to the untreated, incubator control (positive viability control) and the Triton X-100 control (cell death control) remained above 90% viable for all test article exposures and air exposures. (b.) Viability of test article-exposed tissues was expressed relative to the matched air-exposed tissues. Viability of tissues exposed to test articles remained ~100% of the matched air-exposed tissues. Mean viability  $\pm$  s.d. is shown (n=3).

### 6c. Barrier function of EpiAirway™ tissues determined by TEER following exposure to test articles.

Barrier function was assessed by TEER following exposure to test article aerosol, air and control treatments. All tissues had pre-exposure TEER values  $>300 \Omega \cdot \text{cm}^2$ , indicative of an

intact epithelial barrier. Following exposure to Triton, the positive control for cell death, barrier function was completely lost (TEER values  $<1\%$  of pre-exposure value,  $\leq 3.4 \Omega \cdot \text{cm}^2$ ). Untreated, incubator control tissues maintained barrier function with TEER values  $>300 \Omega \cdot \text{cm}^2$ . Likewise, all air-exposed and test article-exposed tissues maintained barrier function (TEER values  $>300 \Omega \cdot \text{cm}^2$ ) after treatment (**Figure 5a**). Barrier function did not differ between tissues exposed to the test article aerosol and tissues exposed to control air, suggesting that the test articles did not impact tissue barrier function (**Figure 5b**).

**Figure 5. Barrier function of control-treated tissues determined by TEER.** Barrier function following treatment was assessed by TEER. (a.) Mean  $\pm$  s.d. absolute TEER values expressed in  $\Omega \cdot \text{cm}^2$  are shown ( $n=3$ ). All tissues exhibited pre-exposure TEER values  $\geq 300 \Omega \cdot \text{cm}^2$ , indicative of an intact epithelial barrier. All post-exposure TEER values were  $\geq 300 \Omega \cdot \text{cm}^2$ , demonstrating that none of the treatments impaired tissue barrier function. (b.) The % pre-exposure TEER values expressed relative to the matched air values show TEER was similar between treatments. Mean  $\pm$  s.d. is shown ( $n=3$ ).


## 6d. Effects of test articles on EpiAirway™ tissue morphology determined by histological assessment.


At the 24-hour time point, all tissues were formalin fixed, paraffin-embedded, sectioned and stained by H&E to assess overall tissue morphology. Untreated, incubator control tissues, air control-exposed tissues and test article-exposed tissues all showed comparable morphology that is typical of a healthy pseudostratified respiratory epithelium (**Figure 6**). No tissue damage or abnormalities were detected. The Triton-treated, negative control tissues were almost entirely destroyed, with only a few cells remaining attached to the insert membrane.


**TA2: Sleep Flavor  
+ Vehicle**


**TA2: Matched Air**


**TA3: Energy  
Inhaler (full  
formulation)**

**TA3: Matched Air**

**Figure 6. Morphology of EpiAirway™ tissues 24 hours post-treatment.** Tissue morphology was assessed by H&E staining of tissue sections (n=3). One representative image from each treatment is shown at a low and high magnification. Images are labeled with the Slide ID number.

## 7. Conclusion

This study evaluated the effects of Nutriair Energy Inhaler and Nutriair Sleep Inhaler, along with the vehicle control aerosol (pure 50:50 PG/VG, Energy vehicle + flavor and Sleep vehicle + flavor) on EpiAirway tissue viability, epithelial barrier function and tissue architecture. Tissue viability, determined using the LDH assay, was unaffected by exposure to the fully formulated test articles and the vehicle control test articles compared to matched air exposure, indicating no cytotoxic effects of the test articles on the EpiAirway tissues at the doses tested. This result was further supported by the evaluation of epithelial tight junction integrity by TEER measurement. TEER showed that all tissues maintained barrier function following exposure to the test articles. Lastly, H&E stained tissue cross-sections demonstrated that the tissue structure and differentiation was similar between untreated tissues and tissues exposed to air, vehicle control aerosol and full formulation aerosol. Follow-up experiments could be conducted to further confirm and extend these results. For example, additional functional endpoints like changes in gene expression or oxidative stress markers could be examined or higher doses could be evaluated to confirm that toxicity would not occur under extreme exposures. Conditioned cell culture media were collected and saved for possible future evaluation of cytokine secretion following test article exposure. In conclusion, the Nutriair Energy Inhaler and Sleep Inhaler did not impair EpiAirway tissue viability, barrier function or morphology under the conditions tested here. EpiAirway™ tissues are a valuable *in vitro* tool for screening and evaluating the effects aerosol products on human respiratory epithelial tissues at the air-liquid interface.

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(3226 Benneth Street North  
St. Petersburg, FL 33713)

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